

ASM


Medicare

ASM
Corner Shower
Stool 02010.
User Instructions.

Picow Farm Road Ind Est
Runcorn. Cheshire.
WA7 4UG.
Tel: 01928 574301.
Fax: 01928 575130.
sales@asmmedicare.com



Corner Shower Stool Code No: 02010.



On receipt of your new **ASM** Corner Shower Stool please check the contents for signs of damage, report any defects to ASM immediately.
The stool comes fully assembled & ready for use and is height adjustable, to adjust please follow the instructions below.
This product is only intended to be used in a shower cubical or similar area to provide a seat whilst showering.

- 1} Lay the stool on it's side and remove the plastic clips from two legs.
- 2} Slide the legs to the required position, ensuring that both legs are adjusted to the same length. Fig 1
- 3} Replace the plastic clips, you may have to twist the leg slightly to allow the steel pin to locate completely, this is **VERY** important & the clips **MUST** be fully engaged.
- 4} Turn the stool over and repeat for the other leg.
- 5} Place the stool back on it's feet, on a level surface and check that it does not rock. If it does readjust the legs until it is stable.
- 6} The stool has a removable seat for ease of cleaning, to remove simply hold the stool & pull the seat free. To replace, it may be necessary to tap the seat with your hand. **Do not** use a heavy object as this may damage your seat.
- 7} To clean the stool a mild proprietary soap and or disinfectant may be used. Dry with a soft clean cloth after cleaning and every use.

Note:

This product should only to be used for it's intended purpose. We do not recommend that the stool is used in a tilting position as this could result in slipping and personal injury may occur.

Maximum User Weight = 190kg / 30st.

For details of all our other products please visit our web site at www.asmmedicare.com

Aug 2005.

Better By Design



IMPORTANT NOTICE

Please read & fully understand before using or allowing others to do so.

