

# ASM

ASM User Instructions.  
Over Bed Tables.

## Medicare

Picow Farm Road Ind Est  
Runcorn. Cheshire.  
WA7 4UG.

Tel: 01928 574301.

Fax: 01928 575130.

e-mail: [sales@asmmedicare.com](mailto:sales@asmmedicare.com)



## ASM Derwent / Essex Over Bed Tables.



On receipt of your new ASM Derwent / Essex Over Bed Table please check the contents for signs of damage, report any defects to ASM immediately.

The tables come fully assembled & ready for use and is height adjustable, to adjust please follow the instructions below. This product is only intended to provide a portable, lightweight rest, suitable for reading, writing and placing meals on. It can be used either over a bed or wheelchair as required.

The Derwent comes with a teak effect laminate table top with either one or two retaining edges. The Essex comes with a ivory coloured laminate table top, again with either one or two retaining edges. Please specify on order.

Both are easy to clean using a damp cloth.

The tables have the facility to adjust the height & angle of the top, to provide a more comfortable reading position, as illustrated.

- 1} To adjust the top release the two handwheels and slide the top to the required height and angle.
- 2} Ensure that both handwheels have been fully re-tightened before placing anything on them.

The tables must not be used as a seat or for personal support, placing your full weight on the table could lead to product failure and possible personal injury.

To clean the tables a mild propriety soap and or disinfectant may be used, dry after cleaning with a soft cloth.

### Note:

This product should only to be used for it's intended purpose.

For details of all our other products please visit our web site at [www.asmmedicare.com](http://www.asmmedicare.com)



October 2009

Better By Design