

# ASM

User  
Instructions  
ASM Wren  
Kitchen /  
Shower Stool  
Heavy Duty  
Range.

## Medicare

Picow Farm Road Ind Estate  
Runcorn. Cheshire. WA7 4UG.  
Tel: 01928 574301  
Fax: 01928 575130  
E-mail: [sales@asmmedicare.com](mailto:sales@asmmedicare.com)

### Wren Kitchen / Shower Stool Heavy Duty Range.



On receipt of your new ASM Wren Stool please check the contents for signs of damage, report any defects to **ASM Medicare Ltd** immediately.

The heavy duty range of Wren Stools come fully assembled and ready for use. They are not height adjustable, but are available in a number options;

Code 02300HD - No Arms or Backrest. {shown}

Code 02302HD - Loop Arms No Backrest.

Code 02303HD - Loop Arms & Backrest.

Code 02304HD - Loop Arms & Padded Backrest. {shown}

- (1) Periodic checks of all screws and fittings is recommended, replace if required.
- (2) To clean the Heron Shower Chair, a mild proprietary soap and / or disinfectant may be used. After cleaning dry with a soft clean cloth.

The Maximum User Weight for these stools is: **190kg / 30st.**

The intended use of this product is to provide the user with a fixed ht “perching” type seat for those unable to sit down fully.

**NOTE:-** These products are only to be used for their intended purpose.  
For details of all our products please visit our web site at [www.asmmedicare.com](http://www.asmmedicare.com)



### IMPORTANT NOTICE

Please read & fully  
understand before  
using or allowing  
others to do so.

